

Product Variants	Product Description	Needle inch	catheter ø outer/inner/length (mm) (G)	Art Code
Complete Sets with LOR Syringe + Filter + Catheter + Needle	Perifix 401	18 x 3 1/4 "	0.45 /0.85 / 1000 (20G)	4514017
	Perifix 402	16 x 3 1/4"	1.05 / 0.60 / 1000 (19 G)	4514025
	Perifix soft 501	18 x 3 1/4 "	0.45 /0.85 / 1000 (20G)	4515013
	Perifix soft 502	16 x 3 1/4"	1.05 / 0.60 / 1000 (19 G)	4515021
Filter Sets with Filter + Catheter + Needle	Perifix 400	18 x 3 1/4 "	0.45 /0.85 / 1000 (20G)	4514009
	Perifix 405	16 x 3 1/4"	1.05 / 0.60 / 1000 (19 G)	4514050
Mini Sets with Catheter + Needle	Perifix 300	18 x 3 1/4 "	0.45 /0.85 / 1000 (20G)	4513002
	Perifix 310	16 x 3 1/4"	1.05 / 0.60 / 1000 (19 G)	4513100
Epidural Catheter	Perifix 315 (18 G)	-	0.45 /0.85 / 1000 (20G)	4513150
	Perifix 325 (16G)	-	1.05 / 0.60 / 1000 (19 G)	4513258
Epidural Needles	Perican G 18	18 x 3 1/4 "	-----	4512383
	Perican G 16	16 x 3 1/4"	-----	4512782
Pediatric Complete sets*	Perifix ONE Paed G 20	20 G x 2"	0.60 x 750 mm (24G)	4512014c

* With LOR Syringe + Filter + Catheter + Needle + Pin Pad

The Aesculap Academy Guidance Center

Our guidance centers provide up to date and comprehensive training programs to doctors in the field of Peripheral Nerve blocks (PNB) & Pain Management using RA techniques.


To download form , you can go to below mentioned link.
<http://www.aesculap-academy.com>

For any further detail you can visit our website - www.aamembersprivy.com
 Or you can contact the following : jagdishchandra.mody@bbraun.com
 contact no.: 022-66682161



Perifix®

On the right track in Epidural Anesthesia



Help improve procedural outcomes and reduce risks

Epidural Anesthesia

Each patient presents a unique anatomical challenge, therefore precision and versatility are equally important in epidural procedures. With B. Braun epidural anaesthesia products you are ensured accurate construction and reliability.

Tuohy-type needles – Perican®

- The unique Tuohy curve and bevel design: Excellent feel and differentiation of tissue layers.
- Highest quality of bevel grinding: Consistent tactile "feedback" in each procedure.
- Polished inner bevel edge: Minimizes the risk of catheter shearing.




Polished inner bevel edge minimizes the risk of catheter shearing



Perican® Paed – Tuohy pediatric epidural needles

Perican® Tuohy epidural needles

Safe catheter insertion using the threading assist guide



User benefits

- Unique high precision Tuohy curve and bevel design for excellent tactile feedback and differentiation of tissue layers
- Superior bending stiffness – perfect control of needle advancement
- Full range of adult and pediatric sizes for all procedures

Perifix® Epidural Kit

User benefits

- Perifix® Standard catheter
 - Excellent shaft stability and kink resistance while still flexible for easy placement and removal
 - Crystal clear for easy visualization of fluids
- Perifix® Soft Tip catheter
 - Additional features to the above:
 - Flexible and atraumatic soft tip design
 - Reduced risk of complications



Perifix® catheter

User benefits

- Perifix® filter
 - Improved patient comfort and hygiene
 - Flexible fixation – no time wasted
- Perifix® connector
 - Safe catheter fixation
 - Easy to use ergonomic design
 - Clear feedback of action "click and ready"



Perifix® epidural flat filter



Perifix® connector



Perifix® LOR syringe

User benefits

- Perifix® LOR syringe
 - Easy-glide plunger with excellent glide characteristics
 - Parabolic graduation markings prevent confusion with conventional syringes
- Latex-free

