

Spinal Anesthesia

Needles for Your Sensitivity

Single-shot technique

For many anesthesiologists, the choice of optimal needle design is a matter of individual preference and tactile sensitivity.

B. Braun therefore offers two different bevel geometries, covering the full range of gauge sizes.

Every needle type has its special characteristics. What unites them is the highest quality and a design optimized for the intended use.

User benefits

- The ideal needle in your experienced hands
- 2 different bevel configuration to match your personal preference and application
- Translucent needle hubs provide rapid detection of CSF flashback

Hub design

The grooved transparent hub provides rapid detection of CSF flashback for confirming accurate placement. Its color-coded stylet ensures rapid identification of needle size and enables perfect matching: each stylet fits exactly into the respective needle's lumen.



Spinocan® with Quincke bevel

- Sharp bevel for low puncture force.
- Safe control due to reliable sensation of the tactile dura click
- Low rate of PDPH thanks to fine gauge needle design.
- Perfect stylet fit, i.e. no tissue coring (risk of epidermoid tumor formation)



Pencan® with pencil point bevel

- Safe handling through easy identification of tissue structures, excellent tactile feedback.
- Protection against PDPH, paresthesias, vascular puncture and nerve damage thanks to the atraumatic tip design.
- Sturdy design, no tip deformation even after bone contact.
- Optimal needle-tip design: Minimal intrathecal penetration depth, avoiding incomplete blocks

The Aesculap Academy Guidance Center



Our guidance centers provide up to date and comprehensive training programs to doctors in the field of Peripheral Nerve blocks (PNB) & Pain Management using RA techniques.

To download form, you can go to below mentioned link.

<http://www.aesculap-academyap.com>


For any further detail you can visit our website - www.aamembersprivy.com


Or you can contact the following : jagdishchandra.mody@bbraun.com


contact no.: 022-66682161

Spinocan®


Pencan®

Spinocan®	Size	Length ϕ (mm)	Needle ϕ (mm)	Length (inch)	Sales unit/Pcs.	Code No. (REF)
	27 G	88	0.42	3½"	25	4503902
	26 G	88	0.47	3½"	25	4502906
	25 G	88	0.53	3½"	25	4505905
	25 G	120	0.53	4¾"	25	4505913
	23 G	88	0.65	3½"	25	4506014
	22 G	88	0.70	3½"	25	4507908
	20 G	88	0.90	3½"	25	4509900
	19 G	88	1.10	3½"	25	4501195
	18 G	88	1.30	3½"	25	4501390

Pencan®	Size	Length ϕ (mm)	Needle ϕ (mm)	Length (inch)	Sales unit/Pcs.	Code No. (REF)
	25 G	88	0.53	3½"	25	4502019
	27 G	88	0.42	3½"	25	4502051

Guide needles	Size	Length ϕ (mm)	Needle ϕ (mm)	Length (inch)
	20 G	35	0.90	1½"

* with guide needle



B. Braun Spinal Anesthesia : The ideal needle in your experienced hands